

COMFORT 6 ZE column



Technical parameters:

External dimensions	222 - 104 cm
Angle of back section	0 - 75
Angle of leg section	0 - 45
Autoregression	back / leg
Foot section	horizontal correct
Height adjustment:	50 - 90 cm
TR/ATR tilt	+15/-15
Safe working load	280 kg



Option:

Bed extension
X-RAY / C-ARMS
Digital scale
LED floor lamps
ACP / LCD
CPR manual

15 cm
ready
ready
ready

ANA-MED



4 PCS CORNER CRASH BUMPERS




FOOT SECTION CORECT TO
HORIZONTAL POSITION



SIDERAIL WITH EMBEDDED
CONTROL SYSTEM



QUICK DETACHABLE BOARD



MANUAL CPR



AUTOREGRESSION



CENTRAL CONTROL SYSTEM



BEST SAFETY SIDERAIL