

BED CONFIGURATION

type	autoregression back / leg	central breake	handle control	ACP control	bed plate 850 mm	bed plate 900 mm	expand 150 mm	X-RAY ready	digital scale	CPR manual
COMFORT 6 ZE column	YES	YES	YES	*LCD	NO	YES	*YES	*YES	*YES	*YES
COMFORT 6 ZE	YES	YES	YES	*YES	NO	YES	*YES	*YES	*YES	*YES
COMFORT 4 ZE	YES	YES	YES	*YES	NO	YES	*YES	*YES	NO	*YES
COMFORT 3 ZE	YES	YES	YES	NO	NO	YES	*YES	NO	NO	*YES
COMFORT 6 ZT	YES	YES	YES	*YES	NO	YES	*YES	NO	NO	*YES
COMFORT 4 ZT	YES	YES	YES	NO	NO	YES	*YES	NO	NO	*YES
COMFORT 3 ZT	YES	YES	YES	NO	NO	YES	*YES	NO	NO	*YES
STANDARD 4 ZE	YES	YES	NO	NO	YES	NO	NO	NO	NO	NO
STANDARD 3 ZE	YES	YES	NO	NO	YES	NO	NO	NO	NO	NO
STANDARD 3 BE	YES	NO	NO	NO	YES	NO	NO	NO	NO	NO
STANDARD 4 ZL	NO	YES	NO	NO	YES	NO	NO	NO	NO	NO
STANDARD 3 ZL	NO	YES	NO	NO	YES	NO	NO	NO	NO	NO
STANDARD 3 BL	NO	NO	NO	NO	YES	NO	NO	NO	NO	NO
STANDARD 1B L	NO	NO	NO	NO	YES	NO	NO	NO	NO	NO
STANDARD 2 BF	NO	NO	NO	NO	NO	NO	NO	NO	NO	NO
COMFORT 3 BL	NO	NO	YES	NO	YES	NO	NO	NO	NO	NO

*option

Decor:



ANA-MED