



# COMFORT 6 ZT light



## Technical parameters:

|                       |                    |
|-----------------------|--------------------|
| External dimensions   | 222 - 104 cm       |
| Angle of back section | 0 - 75             |
| Angle of leg section  | 0 - 45             |
| Autoregression        | back / leg         |
| Foot section          | horizontal correct |
| Height adjustment:    | 39 - 70 cm         |
| TR/ATR tilt           | +15/-15            |
| Safe working load     | 250 kg             |



Plastic bed plate  
CPR manual  
Handset with lock  
Lifter