


# MATTRESS CONFIGURATION



Type	Foam	L	W	H	urineproof	waterproof	silver ions	permeable	reversible	gofer
A 01	polyurethane	200	85	8	yes	yse	option	yes	yes	no
A 02	polyurethane	200	85	12	yes	yes	option	yes	yes	no
A 03	polyurethane	200	90	10	yes	yes	option	yes	yes	no
A 04	polyurethane	200	90	12	yes	yes	option	yes	yes	no
A 05	polyurethane	200	90	15	yes	yes	option	yes	yes	no
A 06	polyurethane	200	90	10	yes	yes	option	yes	no	yes
A 07	polyurethane	200	90	12	yes	yes	option	yes	no	yes
A 08	polyurethane	200	90	15	yes	yes	option	yes	no	yes

